

WoW Cart Maintenance & Safety

A. Cleaning and Maintenance

1. The hospital Infection Prevention Team should be consulted with any queries regarding cleaning procedures and processes.
2. To avoid risk of electric shock, do not expose electrical components to water, cleaning solutions or other potentially corrosive liquids or substances.
3. The mouse and keyboard are waterproof and can be washed. Do not immerse Cart or Cart components in liquid or allow liquids to flow into the Cart.
4. Do not use flammable cleaners on Cart.
5. All paints and plastic components on Cart can be cleaned by most commonly used, diluted, non-abrasive solutions, e.g. quaternary ammonia compounds, ammonia enzyme cleaners, bleach or alcohol solutions.
 - It is necessary to clean the surface of the device using hospital approved disinfecting wipes.
 - Clean mouse, keyboard, handles and work surface between patients.
 - Iodine stains can be removed with commonly used cleaners and a soft cloth.
 - Never use steel brush or other abrasive materials that will damage the surface finish.
 - Do not use strong solvents such as trichloroethylene and acetone. These solvents will damage the surface finish. It is recommended that any cleaning solution be tested on a small inconspicuous area to make sure desk surface being not harmed.
6. Adjustment, Service, Replacement - DO NOT attempt to adjust, service or replace any part of the Cart unless directed to do so through approved documentation released by Original Maker (i.e. installation instructions, etc.). Only Original Maker or a Maker-certified Entity may adjust, service or replace Cart components. If any component on the Cart is missing or damaged, the Cart must not be used. Contact Original Maker or your vendor immediately to request a replacement part and corresponded service.
7. Cables - Keep cables neatly organized on the Cart (a variety of solutions are provided with your cart for this purpose). Excess cables should be routed away from moving components with cable clips or equivalent. Review Cable Routing section of this guide, or contact Maker or your vendor for more information.
8. Casters - Check casters periodically to make sure they are clean and free of debris that would prevent cart from smooth traveling. Avoid moving Cart across uneven, dirty or damaged surfaces.
9. Customer Equipment - Make sure equipment being balanced and mounted securely to Cart. Do not reposition Cart components on Main Lift or Main Pole unless instructed to do so in the installation instructions. Moving Cart components too high or too low on the Main Lift may create an unstable condition, leading to equipment damage or even personal injury. Contact Original Maker or your vendor for information about moving Cart components.

B. Safety Alerts Associated with this Product

The following Warnings/Cautions appear in this reference guide or on the cart:

NOTE: Failure to follow on these guidelines may result in equipment damage or personal injury.

CAUTION: DO NOT load Life module with equipment more than the maximum weight capacity in total specified by Original Maker. Make sure optimum life function by testing and if necessary, re-adjusting tension whenever the weight mounted to the Lift changed (i.e., equipment is removed or added). See "Set Lift Tension" adjustment instructions.

CAUTION: DO NOT operate the Cart with missing or damaged components! Do not remove, modify or substitute Cart components without consulting Original Maker or your vendor. If you encounter problems with Cart on installation or operation, please contact Original Maker or your vendor immediately.

CAUTION: DO NOT overtighten fasteners. Overtightening may cause damage to your equipment.

WARNING: Stored Energy Hazard: The work-surface lift mechanism is under tension and will move up rapidly, on its own, as soon as attached equipment is removed. For this reason, DO NOT remove equipment unless the work surface has been moved to the highest position on the tower!

Failure to follow this instruction may result in serious personal injury and/or equipment damage! When shipping the cart, set the work surface lift mechanism to the lowest tension setting.

CAUTION: DO NOT loosen, tighten or remove any other nuts or bolts on the Lift or top of Main Pole. Tampering with nuts or bolts may result in an unstable Cart, leading to equipment damage and/or personal injury.

CAUTION: Release Lift Brake before moving work-desk surface! Moving work-desk surface while Lift Brake is engaged may cause serious damage to Lift Engine.

WARNING: In the event that repair of the Cart is needed, contact Original Maker immediately. Cart repair can only be performed by Original Maker or by an authorized agent approved by the Original Maker.

WARNING: This cart is not intended for use in a flammable, anaesthetic mixture or oxygen rich environment.

C. Configuration & Safety

1. Do Not Connect Additional multiple socket-outlet or extension cord to the medical system.
2. When used in a Medical Electrical system, connect only equipment that complies with IEC, ISO, UL/ANSI, or CSA standards that are relevant to that equipment.
3. Risk of shock or personal injury when connecting non-medical equipment, supplied as part of a system, directly to the wall outlet when non-medical equipment is intended to be supplied by the multiple socket outlet.
4. Risk of shock or personal injury when connecting any equipment that has not been supplied as part of the medical system, to the multiple socket outlets.