





EMR Quick Reference Guide

Access and Set Up – Customising a Navigator Band for Interactive View and Fluid Balance (iView)



1. Click on **Interactive View and Fluid Balance (iView)** from Table of Contents (TOC). In the **Navigator Bar** at the top of the page select **View > Layout > Navigator Bands**.



2. The **Navigator Documents** window will open, from the panel on the left **select the band** you wish to bring through, in this case the **Adult Sepsis Pathway** and **click on the →**.



- The selected band will now appear in the panel on the right. Click on **OK**
- Changes to Navigator Bands window will appear advising that you will need to close the patient chart so you can view the changes. Click on **OK**.



- Once you have reopened the patient chart go to **iView** and the band you added will appear