

Orders - Intravenous Cannula Care Task



This Quick Reference Guide will explain how to: Order - Intravenous Cannula Care Task

- This will generate an order each shift to document cannula care, and will automatically cease after 72 hours.

- Go to **Menu** → **Orders** → **Add** → Search “Intravenous Cannula Care”

Search: intravenous | Advanced Options | Type: Inpatient

- Intravenous Cannula Care
- Intravenous Pyelogram CT
- Nu... Cholangiogram Intravenous CT

US Doppler Abdomen

- Fill in the details as - **Frequency: TDS, Duration: 72, Duration unit: hr(s).**

This will generate an order each shift to document cannula care, and will automatically cease after 72 hours.

Details for **Intravenous Cannula Care**

Details | Order Comments | Diagnoses

*Requested start date and time: 19/08/2019 1105 AEST | Frequency: TDS

Duration: 72 | Duration unit: hr(s)

Stop date and time: 22/08/2019 1359 AEST | Constant order: Yes No

PRN: Yes No | Special instructions:

- It's recommended if you are placing this order regularly to save it to your favourites, at this point. Please see the **QRG Access and Set Up – Creating a Favourites Folder** for the full favourite's workflow.
- Click **Sign** in the bottom right corner.

Sign