Continuous Infusions

Continuous Infusions = both "truly continuous" and "bag by bag" infusions.

They display under the **Continuous Infusions** section of the MAR/MAR Summary/ Orders page:

Continuous Infusions

amiodarone (additive) 300 mg
Glucose (PVC free) 5% infusion (BAG BY BAG) 100 mL
100 mL, IV Infusion, Rate: 200 mL/hr, Duration: 1 bag(s), First dose
07/09/2020 12:42:00, Indication: Loading dose, Stop date 07/09/2020
13:11:00, Use PVC free bag., Total volume (mL): 100
Refer to WH procedure. Given after loading dose

Only the RATE of a continuous infusion will display on the MAR/MAR Summary.

Note: The infuse over time is automatically calculated based on the rate and the bag volume. In some situations, this may mean the infuse over time is clinically inappropriate.



