## Infusions- Pausing an Infusion -Nursing/Midwifery



Digital Health Quick Reference Guide

## Key message

- If the infusion is required to be temporarily paused, manually stop the IV infusion pump using the Rate Change function.
- 1. Navigate to the MAR and click on the Pending tile to open up the MAR Charting Window.
- 2. Click on **Rate Change** and **modify** the rate in the bottom section to **0ml/hr** to indicate the infusion is paused.
- 3. Click Apply
- 4. Discern alert will display indicating the rate change differs from the prescribed order.
- 5. Select Yes.
- 6. Add a comment if required in MAR Charting window
- 7. To Sign off the MAR Charting window click on the Green tick
- 8. Refresh the MAR and review the infusion rate change displays correctly

**NB:** The End bag tile appears at the expected time the infusion bag is empty (based on the rate and documented volumes infused). Given the rate is 0ml/hr the End Bag will appear on the right side of the MAR pending tile.

9. When infusion can be recommenced repeat steps
a) – c)

Time View	Medications	21/10/2020 02:14	21/10/2020 02:13	21/10/2020 02:12
Scheduled	sodium chloride 0.9% infusion (BAG BY BAG) 1.000 mL	Pending ast bag started:	End Bag	
Unscheduled	1,000 mL, IV Infusion, Rate: 125 mL/hr, Duration: 2 bag(s), First	1/10/2020 02:11 AEDT	1 of 2	
PRN	volume (mL): 1,000			
Continuous Infusions	Administration Information		* Rate Change 0 mL/hr	25 mL Infuse
	Sodium Chloride 0.9%		*	125 mL



